




11 sphered jade projector



DESCRIPTION OF THE CONTROL PANEL FEATURES:

1. Plug in the projector into the outlet.
2. Turn on the projector by pressing on power button.
3. Buttons to control temperature are called temperature add-to increase temperature, and temperature cut- to decrease the temperature. You can control the needed temperature by clicking on the add and cut buttons. The temperature ranges from 20 to 70 degrees Celcius. The temperature goes down by 1 degree.
4. Button time-by pressing on this button the time can be controlled. The time for the procedure ranges from 10 minutes to 1 hour and 20 minutes. The vibration mode lasts for 10 minutes, it will turn off automatically after 10 minutes.
5. Vibrate button: turn it on and enjoy the vibration mode for 10 minutes at a time.
6. The negative ion button-turn on the button to start getting negatively charged ions in the air.

This section contains safety instructions to help you avoid the risk of body injuries.

Not adhering to these rules makes the warranty not applicable!

1. Use the product on a flat, hard surface.
2. Disconnect the cord from the outlet if you are not using the product any longer.
3. Do not bend the cord or place heavy objects on it.
4. Do not cover the control panel with a blanket when using the device.
5. Check to see if the remote control panel is connected properly to the device before plugging the cord into the power outlet.
6. Don't use the easy flammable products such as candles, don't smoke near the device while it's plugged in.
7. Don't use needles or any other metallic devices near the product, don't let the metallic devices get in contact with the product.
8. Don't use the product in conjunction with other devices that heat up.
9. Don't allow children to use the product without adult supervision.
10. Don't attempt to open up the product and fix it by yourself.
11. If you suspect any issues with the product, stop using it and call the professional to get it fixed.
12. Don't use any chemicals to clean the product. Only use a soft dry piece of material to clean it.
- 13. If you have a pacemaker or any kind of implant in your body, make sure to consult a doctor prior to starting to use the product.**

Calcium, magnesium and iron silicates with very small amphibol pieces are included in the jade.

Since ancient times, the jade stone was considered the first remedy for all diseases. People wore it on their necks as a treatment and to prevent any diseases to start. In various medical books published in the East, Jade is described like a stone that heals damaged organs and affects the tissues of the body for faster healing. In China, many diseases are successfully treated with jade powder.

Digestive system diseases are treated with jade stones. Pieces of this mineral were put on sore spots in person's body to reduce pain after injuries, any deformities, or arthritis.

In many countries of the world jade balls are now used to massage the body and face. Jade has the ability to smooth the skin of the face, so beauticians, finishing the massage, put pieces of jade on the faces of women. Among other things, stone is an excellent remedy for headaches, dizziness, atherosclerosis, strokes and eye diseases. It eliminates toothache, renal colics, normalizes cardiovascular activity, improves the work of the urinary system. Perfectly copes with nervous diseases, has a calming property, helps to get out of depression. The stone helps to restore the circulatory system and bone marrow.

Ionization

While vacationing at the sea, walking through the forest or being in the mountains, we feel how our bodies become enriched with oxygen.

This is because negatively charged particles have the ability to quickly enter the body and have a beneficial effect on blood circulation.

Nature is filled with negatively charged ions: in the forest 200 - 1000 ions/cm³, on the sea 2000 ions/cm³, in the mountains 10,000 ions/cm³, at a waterfall of 18,000 ions/cm³.

When jade is heated in the projector, the air becomes saturated with negatively charged ions. The presence of negatively charged ions in the human body is necessary, as they are inevitable participants in metabolism, improve the permeability of cell membranes, help to burn fat cells. From muscle contractions, conductivity excite pulses along nerve fibers. Ionization also has a beneficial effect on the maintenance of neural self-regulation, helps to activate water molecules in the body by normalizing the viscosity of the blood.

Reducing the amount of negatively charged ions by significantly worsening the process of cell absorption of nutrients and removing toxins from the body. Only with the balance of two types of ions (positive and negative) can the norm be maintained by vital forces and support of the body with infections.

Negatively charged ions enter the body when a person uses the equipment. It

использует оборудование

Provides the effect of:

1. Anti-allergic action;
2. Normalization of blood pressure;
3. Improving cardiac activity;
4. Reducing spasms in the musculature of lungs, slowing down and increasing breathing patterns.
5. Improving the functions of the mucous membranes and airways;
6. Strengthening immunity;
 1. Destruction of free radicals (which contribute to the onset of cancer);
 2. Normalization (red blood cell accession ratio);
 3. It has a soothing effect on the nervous system.

But it does not happen suddenly, we have to arm ourselves with patience and constantly, and use the product for a few months to be able to see the results.

“Universal doctor.”

11-ball projector is made like this: Ten 10 balls are arranged equally and the eleventh is concave. This allows you to place it on biologically active points on the back, but without touching the spine of the back. **11 jade projector** is versatile to use.

You can use the jade projector to massage your children, you can use it on their necks and ears and feet.

It is no exaggeration to say that the compact **Jade Projector** is able to cope with many, even very serious problems. Joint diseases: It is recommended to use the projector with the heating action to relax the muscles prior to going to bed. When you heat up the muscles spasms decrease. When you heat up the joints of the body they start healing because they get enough oxygen which is a food for the joints to be able to heal. It is recommended to use projector on your joints in the morning and at night.

Getting rid of the toxins in the body.

Recent studies have shown that infrared diffusion stimulates the removal of poisonous substances, including lead and mercury, from the body's cells, they come out with sweat and urine.

Cleansing the body of toxins is a prerequisite for preventing various diseases and increasing good health.

Pain decreasing effect

Infrared heat helps to fight pain in muscles and joints, eliminating cramps and stiffness in the limbs.

Anti-inflammatory effect

The jade projector will help effectively eliminate inflammation of viruses of any localization. During cold, the reproduction of pathogenic bacteria and viruses in the body is suppressed. The use of the device will significantly increase the recovery time or avoid colds.

Jade is simply irreplaceable in inflammation of ear, nose and throat and with respiratory organs. The use of equipment calms down the nervous system, helps to eliminate insomnia, decreases stress and anxiety.

The use of jade makes it possible to get a wonderful cosmetic effect. Under the influence of infrared rays, blood circulation in the skin improves. The pores of the skin are cleaned, dead cells are removed, the skin becomes smooth and elastic. Scars on the skin soften, and in some cases dissolve. The whole complexion looks better.

Hypertension and atherosclerosis

Warming and massaging the neck and shoulder girdles with an 11-ball projector decreases the spasms of the cervical muscles and restores blood flow in the spinal arteries.

To improve the blood supply to the head, application of massager to the carotid arteries is needed, the loss of capillaries is restored, and white blood cells "eat" atherosclerotic plaques.

It is recommended to use projector to the stomach area-in the liver region, in the bottom of the stomach for 15-30 minutes with the temperature 40-55 degrees Celcius, it can be moved from one region to another depending on the area of where the pain is.

The projector can be applied to the shoulder region or any other joint in the body for 15-30 minutes at a time depending on the pain localization.

Massage of the face and neck with the use of jade projector can be done at 40-50 degrees.

Massaging the muscles of the back, legs, arms and buttocks at 55-65 degrees.

Manual projector can be used to heat up biologically active points in the body with the temperature starting at 55 degrees for two minutes at one point.

To decrease and prevent burns it is recommended to hold the massager at one area for no longer than 15 minutes, do not place the massager on the bare skin.

When you have a headache it is recommended to place the massager on the top of the head for 2-3 minutes, it can be placed on the temples and on forehead.

Ringling in the ears: it is recommended to fold the ear in and hold the projector next to it for 2-3 minutes.

Runny nose: on top of nose and in between the eyebrows.

Sore throat, dry cough-next to the chin.

Teeth pain: underneath the jaw and where the pain is.

Bronchitis, licer, stomach, large intestine- next to the organs where the pain is.

Constipation, inability to get pregnant, urinary tract problems-the bottom of the abdomen.

Knee joint pain- on the side of the knee and under the knee.

Caution!

Those who are predisposed to increased blood pressure are advised to start sessions at 40 degrees.

At a "moderate" temperature of 40-45 degrees jade can be used by young children, people suffering from heart disease, kidney disease, urinary system, stroke or heart attack.

Those who suffer from chronic problems should start their therapy with a low temperature. Each week, for example, the temperature can be raised by 5 degrees. Over time, the degrees may increase. The higher the temperature, the stronger the healing effect will be.

Reaction of the body to the recovery process

Healing while using the equipment cannot be done without aggravation of the symptoms. But it should not bother you, it is very important that adverse reactions to the equipment are treated with understanding, because it is partly associated with the transition from a chronic form of a disease or a problem to an acute form. The reaction to recovery is not only felt in perfectly healthy people.

How to react to adverse reactions:

1. In the period of aggravation of the symptoms, jade equipment can still be used (low temperature procedures can be applied).
1. Don't forget to drink clean water (non-carbonated) for an average of 2 litres, especially before breakfast and before meals/during the day to help the body cleanse.
1. Eat more fruits and vegetables, limit fast food, fatty and spicy foods.